

**X-RAY
ACCESSORIES
DARKROOM LAMP**

DARKROOM LAMP *DL1*

***Can be used on both
regular and ortho x-ray film.
Can be installed with a single screw***



■ Features


Long service life : Troublesome lamp replacements are not required with LED type.

High brightness : With LED type, darkroom work can be easily preformed for a long period.

Easy installation : The lightweight design allows installation with a single screw.

Common type : Can be used for both regular and ortho x-ray film.

■ Dimensions



■ Specifications

Type	DL1
System	Special high-brightness LED
Protection	Waterproof (Moisture-proof type)
Size	120 (Dia.) × 160 (H) (Main part 94) mm
Weight	300g
Power	AC220/240 2W

⚠ Safety precautions • Prior to use, please make sure to read the "Precautions for use"